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BIBLIOGRAPHIE

- Free, A.H and Free, H.M.: Urinalysis, Critical Discipline of Clinical Science. CRC Crit. Rev. Clin. Lab. Sci. 3(4): 481-531; (1972).
- Yoder, J.Adams, E.C., and Free. H.M.: Simultaneous Screening for Urinary Occult Blood, Protein, Glucose and pH. Amer. J. Med Tech. 31:285; (1965).
- Tietz, N.W.: Clinical Guide to Laboratory Tests; W.B. Saunders Company, (1976).
- Burtis, C.A. and Ashwood, E.R.: Tietz Textbook of Clinical Chemistry 2nd Ed. 2205; (1994).
- Shcherstin, B. and Fritz, H.: Subnormal Levels of Glucose in Urine. JAMA 201:129-132; (1967).
- McGarry, J.D.: Lilly Lecture, 1978: New Perspectives in the Regulation of Ketogenesis. Diabetes 28: 517-523 May, (1978).
- Williamson, D.H.: Physiological ketoses, or Why Ketone Bodies? Postgrad. Med. J. (June Suppl.): 371-375: (1971).
- Paterson, P. et al.: Maternal and Fetal Ketone Concentrations in Plasma and Urine. Lancet: 862-865; April 22, (1967).
- Fraser, J. et al.: Studies with a Simplified Nitroprusside Test for Ketone Bodies in Urine, Serum, Plasma and Milk. Clin. Chem. Acta II: 372-378; (1965).
- Henry, J.B. et al.: Clinical Diagnosis and Management by Laboratory Methods, 16th Ed. Philadelphia: Saunders; (1979).
- Hayashi A et al:Journal of Japan Diabetes Society, 35,819,1992

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Yercon Diagnostic Co., Ltd.

No.2 Building, Ledong Industrial

Estate 3248 Century Ave. Tech. & Economic Development Zone

Changchun P.R. China 130032

E-Mail: yerconlab@yerconlab.com

Tel: +

REAGENT STRIPS FOR URINALYSIS

(For both professional & self-testing use)

This package insert is to be used with Reagent Strips for Urinalysis (from 1 parameter to 11 parameters). For the semi-quantitative and qualitative detection of Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocytes and Ascorbic Acid in urine. For self-testing, in vitro diagnostic use. Single use only; do not reuse.

SUMMARY

Reagent strips for urinalysis (URS) are firm plastic strips to which several different reagent areas are affixed. Depending on the product being used, reagent strips for urinalysis provide tests for Glucose, Bilirubin, Ketone (Acetoacetic acid), Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocytes, and Ascorbic Acid in urine. Test results provide information regarding the status of carbohydrate metabolism, kidney and liver function, acid-base balance, and bacteriuria. Please refer to the color chart for the specific test parameters of the product you are using.

Reagent strips for urinalysis are packaged along with a drying agent in a pouch with a twist-off cap. Each strip is stable and ready to use upon removal from the pouch. The entire reagent strip is disposable. Results are obtained by direct comparison of the test strip with the color blocks printed on the color chart. No calculations or laboratory instruments are required.

TEST PRINCIPLE

Glucose: This test is based on a double sequential enzyme reaction. One enzyme, glucose oxidase, catalyzes the formation of gluconic acid and hydrogen peroxide from the oxidation of glucose. A second enzyme, peroxidase, catalyzes the reaction of hydrogen peroxide with potassium iodide chromogen to oxidize the chromogen to colors ranging from yellow to greenish-blue and dark blue.

Bilirubin: This test is based on the coupling of bilirubin with a diazotized dichloroaniline in a strongly acid medium. The colors range from light tan to reddish-brown.

Ketone: This test is based on the reaction of acetoacetic acid with sodium nitroprusside in a strongly basic medium. The colors range from beige or buff-pink color for a "Negative" reading to pink and pink-purple for a "Positive" reading.

Specific Gravity: This test is based on the apparent pKa change of certain pretreated polyelectrolytes in relation to the ionic concentration. In the presence of an indicator, the colors range from dark blue or blue-green in urine of low ionic concentration to green and yellow-green in urine of higher ionic concentration.

Blood: This test is based on the pseudoperoxidase action of hemoglobin and erythrocytes which catalyze the reaction of 3,3', 5,5'-tetramethyl-benzidine and buffered organic peroxide. The resulting colors range from orange to yellow-green and dark green. Very high blood concentration may cause the color development to continue to dark blue.

pH: This test is based on the well known double pH indicator method, where bromothymol blue and methyl red give distinguishable colors over the pH range of 5.0-9.0. The colors range from red-orange to yellow and yellow-green to blue-green.

Protein: This test is based on the protein error-of-indicator principle. At a constant pH, the development of any green color is due to the presence of protein. Colors range from yellow for a "Negative" reaction to yellow-green and green to blue-green for a "Positive" reaction.

Urobilinogen: This test is based on a modified Ehrlich reaction in which p-diethylaminobenzaldehyde reacts with urobilinogen in a strongly acid medium. Colors range from light pink to bright magenta.

Nitrite: This test depends on the conversion of nitrate to nitrite by the action of Gram-negative bacteria in the urine. The nitrite reacts with p-arsanilic acid to a diazonium compound in an acid medium. The diazonium compound in turn couples with 1,2,3,4-tetrahydrobenzo(h) quinoline to produce a pink color.

Leukocytes: This test is based on the action of esterase present in leukocytes, which catalyzes the hydrolysis of an indoxyl ester derivative. The indoxyl ester liberated reacts with a diazonium salt to produce a beige-pink to purple color.

Ascorbic Acid: The composition comprises of a complex chelating agent with a polyvalent metal ion in its higher state and an indicator dye that reacts with the metal ion in its lower state to produce a color change from blue-green to yellow.

REAGENTS (Based on dried weight at time of impregnation)

Glucose: 16.3% w/w glucose oxidase (*Aspergillus niger*, 1.3IU); 0.6% w/w peroxidase (horseradish, 3300 IU); 7.0% w/w potassium iodide; 76.1% w/w buffer and non-reactive ingredients.

Bilirubin: 0.4% w/w 2,4-dichloroaniline diazonium salt, balanced with buffer and non-reactive ingredients.

Ketone: 7.7% w/w sodium nitroprusside balanced with buffer and non-reactive ingredients.
Specific Gravity: 2.8% w/w bromothymol blue, 69.0%; poly(methyl vinyl ether/maleic anhydride); 28.2% sodium hydroxide

Blood: 6.6% w/w cumene hydroperoxide; 4.0% w/w 3, 3', 5, 5'-tetramethylbenzidine; 89.4% w/w buffer and nonreactive ingredients.

pH: 0.2% w/w methyl red; 2.8% w/w bromothymol blue; 97% w/w nonreactive ingredients.

Protein: 0.3% w/w tetrabromophenol blue; 99.7% w/w buffer and nonreactive ingredients.

Urobilinogen: 2.9% w/w p-diethylaminobenzaldehyde balanced with buffer and nonreactive ingredients.

Nitrite: 1.4% w/w p-arsanilic acid, balanced with buffer and nonreactive ingredients.

Leukocytes: 0.4% w/w indoxyl ester derivative; 0.2% w/w diazonium salt; 99.4% w/w buffer and nonreactive ingredients.

Ascorbic Acid: 5.8% w/w ferric chloride; 4.9% w/w DTPA; 1.2% dipyriddy; 89.1% w/w buffer and nonreactive ingredients.

WARNINGS AND PRECAUTIONS

- Reagent Strips for Urinalysis are for in vitro diagnostic use. Do not touch the test areas of the Urine Reagent Strips.
- Package contents may be disposed of in the normal household waste after use.
- All package contents are non-toxic and safe when used as directed.
- Do not use if the product seal is broken.
- As with all diagnostic tests, a definitive diagnosis should not be based on the result of a single test, but should only be made by your physician after all clinical and laboratory findings have been evaluated.

STORAGE

- Store at room temperature between 2–30°C (36–86°F) and out of direct sunlight.
- Do not use after the expiration date printed on the label.
- The expiration date is printed on the packaging and the color chart.

SPECIMEN COLLECTION AND PREPARATION

Collect urine in a clean container and test as soon as possible. The use of urine preservatives is not recommended. If testing cannot be performed within one hour after voiding, refrigerate the specimen immediately. Allow a refrigerated specimen to return to room temperature before testing.

MATERIALS

Material provided


Each box contains

- One PET Bottle / Individual pouch containing test strips
- Desiccant (in PE bottle or pouch)
- One color chart (affixed to the bottle, pouch, or box)


Material required but not provided

- Timer
- Sample container
- Disposable glove


TEST PROCEDURE


- Remove one strip from the bottle or pouch and use it immediately.


- Completely immerse the reagent areas of the strip in fresh urine, then remove it immediately to prevent the reagents from dissolving.

- While removing the strip, touch its edge against the rim of the urine container to remove excess urine.

- Blot the long edge of the strip on an absorbent paper towel to remove any remaining urine and prevent cross-contamination between reagent pads.

- Compare each reagent area with its corresponding color block on the color chart and read the results at the specified time. Accurate timing is essential for reliable results.

RESULTS

Results are obtained by direct comparison of the color blocks on the test strip with the color chart. The color blocks represent nominal values; the actual values will vary around the nominal values.

LIMITATIONS OF THE PROCEDURE

As with all tests, definitive diagnostic or therapeutic decisions should not be based on any single test result or method.

Glucose: Moderate amounts of ketone bodies (40mg/dL or greater) may decrease color development in urine containing small amounts of glucose (75–125 mg/dL). However, such a concentration of ketone simultaneously with such a glucose concentration is metabolically improbable in screening. The reactivity of the glucose test decreases as the SG and/or ascorbic acid of the urine increases. Reactivity may also vary with temperature.

Bilirubin: Reactions may occur with urine containing large doses of chlorpromazine or rifampin that might be mistaken for positive bilirubin. Indican (indoxyl sulfate) and metabolites of Lodine may cause false positive or atypical color; ascorbic acid (25mg/dL or greater) may cause false negative results.
Ketone: Color reaction that could be interpreted as "positive" may be obtained with urine specimens containing MESNA or large amounts of phenylketones or L-dopa metabolites.

Specific Gravity: The chemical nature of the specific gravity test may cause slightly different results from those obtained with the specific gravity methods when elevated amounts of certain urine constituents are present. Highly buffered alkaline urine may cause low readings relative to other methods. Elevated specific gravity readings may be obtained in the presence of moderate quantities (100–750 mg/dl) of protein.

Blood: The sensitivity of the blood test is reduced in urine with high specific gravity and/or high ascorbic acid content. Microbial peroxidase, associated with urinary tract infection may cause false positive reactions.

pH: If proper procedure is not followed and excess urine remains on the strip, a phenomenon known as "running over" may occur, in which the acid buffer from the protein reagent area runs onto the pH area,

causing a false lowering in the pH result.

Protein: False positive results may be obtained with highly alkaline urine. Contamination of the urine specimen with quaternary ammonium compounds may also produce false positive results.

Urobilinogen: The test area will react with interfering substances known to react with Ehrlich's reagent, such as porphobilinogen and p-aminosalicylic acid. This test is not a reliable method for the detection of porphobilinogen. Drugs containing azo-dyes (e.g. Azo Gantrisin) may give a masking golden color. The absence of urobilinogen cannot be determined with this test.

Nitrite: The pink color is not quantitative in relation to the number of bacteria present. Any degree of pink coloration should be interpreted as a positive nitrite test suggestive of 10 or more organisms/mL. There are occasional urinary tract infections from organisms, which do not contain reductase to convert nitrate to nitrite.

Leukocytes: Highly colored urine and the presence of the drugs cephalixin (Keflex) and gentamicin have been found to interfere with this test. High urinary protein of 500 mg/dl or above diminishes the intensity of the reaction color. Elevated glucose concentration or high specific gravity may cause decreased results.

Ascorbic Acid: Strong reduction of the metal particles, such as ferrous, stannous, cuprous sulfite, etc. will determine the value high.

EXPECTED VALUES

Glucose: Small amounts of glucose are normally excreted by the kidney. Concentrations as little as 0.1 g/dl glucose, read either at 10 or 30 seconds, may be significantly abnormal if found consistently. At 10 seconds, results should be interpreted qualitatively, for semi-quantitative results, read at 30 seconds only. This test cannot replace patients periodically performing formal medical inspection or diagnostic in a hospital or lab. This test only has the function of routine daily urine glucose measurement.

If trace amounts are found or the urine is positive, it suggests potential clinical conditions such as renal tubular illness, pancreatitis, or diabetes mellitus. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Bilirubin: Normally, no bilirubin is detectable in urine by even the most sensitive method. Even trace amounts of bilirubin are sufficiently abnormal to require further investigation. Atypical colors (colors produced which are different than the negative or positive color blocks shown on the Color Chart) may indicate that bilirubin derived bile pigments are present in the urine sample and are possibly masking the bilirubin reaction.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as obstructive jaundice, chronic liver cirrhosis or common bile duct stones. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Ketone: Normally, no ketones are present in urine. Detectable levels of ketone may occur in urine during physiological stress conditions such as fasting, pregnancy, and frequent strenuous exercise. In starvation diets, or in other abnormal carbohydrate metabolism situations, ketones appear in the urine in excessively large amounts before serum ketones are elevated.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as diabetic acidosis, starvation, or hyperthyroidism. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Specific Gravity: Random urine may vary in specific gravity from 1.003-1.040+. Twenty-four hour urine from normal adults with normal diets and normal fluid intake will have a specific gravity of 0.016-1.022 in severe renal damage the specific gravity is fixed at 1.010, the value of the glomerular filtrate.

If high-density urine is detected, it may indicate clinical conditions such as dehydration or renal insufficiency, possibly even renal adenocarcinoma. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Blood: Any green spots or green color developing on the reagent area within 40 seconds is significant and the urine should be examined further. Note: For females, undergoing this test during menstruation may yield false-positive results. Therefore, we advise against conducting the test during menstruation.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as colonic diverticulum, prothrombin damage, or renal calculus. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

pH: newborn: 5-7 thereafter: 4.5-8 average: 6.3

Protein: In 24-hour urine, 1-14 mg/dl of protein may be excreted by the normal kidney. A color matching any color block greater than trace indicates significant proteinuria. For urine with high specific gravity, the test area may most closely match the trace color block even though only normal concentrations of protein are present. Clinical judgment is needed to evaluate the significance of trace results.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as nephritic syndrome or pyelonephritis. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Urobilinogen: In a healthy population, the normal urine urobilinogen range obtained with this test is 0.2-1.0 Ehrlich Unit/dL. A result of 2.0 EU/dl may be of clinical significance and the same patient sample should be evaluated further.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as bile duct obstruction, intestinal flora imbalance, or extravascular hemolysis. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Nitrite: Normally no detectable amount of nitrite is present in urine. The nitrite area will be positive in a proportion of cases of significant infection, depending on how long the urine specimens were retained in the bladder prior to collection. Retrieval of positive cases with the nitrite test range from as low as 40%, in instances where little bladder incubation occurred, to as high as 80% in instances where a minimum of 4 hours incubation occurred.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as bacteriuria. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Leukocytes: Normal urine specimens generally yield negative results with this test. A trace result may be of questionable clinical significance and it is recommended that the test be repeated using a fresh sample from the same patient. Repeated trace and positive results are of clinical significance.

If a positive result is obtained in the urine test, it suggests potential clinical conditions such as renal calculus, nephritis, urethritis, or cystitis. It's advised to consult a doctor for a definitive diagnosis at a medical facility.

Ascorbic Acid: The daily urinary output of ascorbic acid varies with the intake: it is approximately half of the intake. The average urinary output ranges from 20-30 mg/day. If ascorbic acid is detected in the urine, stop taking ascorbic acid for 24 hours and retest.

False negative and weak reaction of glucose, blood and bilirubin may be observed if:

- Glucose: more than 50 mg/dl ascorbic acid in the sample.
- Bilirubin: more than 50 mg/dl ascorbic acid in the sample.
- Blood: more than 10 mg/dl ascorbic acid in the sample.

Expected Concentration Range

The typical concentration ranges and linear ranges for each parameter of the urine test strips:

Test Items	Test Unit	Expected Test Concentrations					
		-	±	1+	2+	3+	4+
Leucocytes	cells/µL	neg.		ca. 10-25	ca. 75	ca.500	
Nitrite		neg.		pos.			
Urobilinogen	mg/dL	norm		2	4	8	12
Protein	mg/dL	neg.	trace	30	100	500	

			5, 6, 6.5, 7, 8, 9				
pH	Blood	cells/µL	neg.	ca. 5-10	ca. 50	ca. 300	
Ascorbic Acid	mg/dL	neg.		20	40		
Spec. Gravity				1.000, 1.005, 1.010, 1.020, 1.025, 1.030			
Ketones		neg.	trace	15	50	150	
Bilirubin		neg.		1	3	6	
Glucose	mg/dL	norm		50	150	500	1,000

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance characteristics of the Reagent Strips for Urinalysis (URS) have been determined both in the laboratory and in clinical tests. Parameters of importance to the user are sensitivity, specificity, accuracy, and precision. Generally, Urine Reagent Strips (URS) have been developed to be specific for the constituent to be measured with the exception of interferences listed above. (See LIMITATIONS OF PROCEDURE)

The accuracy of visually read strips depends on how well the color blocks on the chart are determined and how effectively the human eye can discern them. Assessing precision in such tests is challenging due to the inherent variability of human vision.

Glucose: This test is specific for glucose; no substances excreted in urine other than glucose are known to give a positive result. The reagent area does not react with lactose, galactose, fructose, or reducing metabolites of drugs; e.g., salicylates and nalidixic acid. This test may be used to determine whether the reducing substances found in urine are glucose.

Bilirubin: This test is considered specific for bilirubin in urine.

Ketone: The ketone test area provides semi-quantitative results and reacts with acetoacetic acid in urine. This test does not react with beta-hydroxybutyric acid or acetone.

Specific Gravity: The specific gravity test permits the determination of urine specific gravity between 1.000 and 1.030. In general, the specific gravity test correlates within 0.005 with values obtained with the refractive index method.

Blood: This test is slightly more sensitive to free hemoglobin and myoglobin than to intact erythrocytes.
pH: The pH test area permits the quantitative differentiation of pH values to one unit within the range of 5-9. The pH reading is not affected by variation in the urinary buffer concentration.

Protein: The test area is more sensitive to albumin than to globulin, hemoglobin, Bence-Jones proteins, and mucoprotein; a negative result does not rule out the presence of these other proteins.

Urobilinogen: The absence of urobilinogen in the urine sample cannot be determined with this test.

Nitrite: Comparison of the reacted reagent area on a white background may aid in the detection of low levels of nitrite ion, which may otherwise be missed. This test is specific for nitrite and will not react with substances normally excreted in the urine.

Leukocytes: This test will not react with erythrocytes or bacteria common in urine.

Ascorbic Acid: This test can detect ascorbic acid in concentrations as low as 0.6 mmol/L in urine.

BIBLIOGRAPHY

- Free, A.H and Free, H.M.: Urinalysis, Critical Discipline of Clinical Science. CRC Crit. Rev. Clin. Lab. Sci. 3(4): 481-531; (1972).
- Yoder, J, Adams, E.C., and Free. H.M.: Simultaneous Screening for Urinary Occult Blood, Protein, Glucose and pH. Amer. J. Med. Tech. 31:285; (1965).
- Tietz, N.W.: Clinical Guide to Laboratory Tests; W.B. Saunders Company, (1976).
- Burtis, C.A. and Ashwood, E.R: Tietz Textbook of Clinical Chemistry 2nd Ed. 2205; (1994).
- Shchersten, B. and Fritz, H.: Subnormal Levels of Glucose in Urine. JAMA 201:129-132; (1967).
- McGarry, J.D.: Lilly Lecture, 1978: New Perspectives in the Regulation of Ketogenesis. Diabetes 28: 517-523 May, (1978).
- Williamson, D.H.: Physiological ketoses, or Why Ketone Bodies? Postgrad. Med. J. (June Suppl.): 371-375; (1971).
- Paterson, P. et al.: Maternal and Fetal Ketone Concentrations in Plasma and Urine. Lancet: 862-865; April 22, (1967).
- Fraser, J. et al.: Studies with a Simplified Nitroprusside Test for Ketone Bodies in Urine, Serum, Plasma and Milk. Clin. Chem. Acta II: 372-378; (1965).
- Henry, J.B. et al.: Clinical Diagnosis and Management by Laboratory Methods, 16th Ed. Philadelphia: Saunders; (1979).
- Hayashi A. et al. Journal of Japan Diabetes Society, 35,819,1992

SYMBOLS

 For in vitro diagnostic use

 Use By

 Consult instruction for use

 Lot number

 Temperature limit for storage

 Do not reuse

 Contains sufficient for <n> tests

Distributed by

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Altonaer Straße 27

10555 Berlin, Germany

Tel: +49 6401 9699 390

Email: service@luxmediq.de

Yercon Diagnostic Co., Ltd.

No.2 Building, Ledong Industrial

Estate 3248 Century Ave. Tech. & Economic Development Zone

Changchun P.R.China 130032

Email: yerconlab@yerconlab.com

Tel: +86 431 84698368 Fax: +86 431 84619986

EC REP

CMC Medical Devices & Drugs S.L.

C/ Horacio Lengo No 18, CP

29006, Málaga, Spain

Tel: +34 951 214 054

Email: sales1@cmcmmedicaldevices.com

YEC-CE-01-010 Rev. 1, 2019-12-09