## One Step Pregnancy Test (Urine)

For self-testing and In Vitro diagnostic only

### INTENDED USE

The LuxmedIQ Pregnancy Urine Test for self-testing is a rapid, one-step lateral flow immunoassay for the qualitative determination of human chorionic gonadotropin (hCG) in urine to aid in the early detection of pregnancy.

It is intended for the over-the-counter use.

### SUMMARY AND EXPLANATION

hCG is a hormone produced by your placenta when you are pregnant. It appears shortly after the embryo attaches to the wall of the uterus. If you are pregnant, this hormone increases very rapidly.

The concentration of hCG increases to 5-50 IU/L one-week post implantation, reaches about 100 IU/L at the time of the first missed menstrual period, and peaks at 100,000-200,000 IU/L at the end of the first trimester. The appearance of hCG soon after conception and its subsequent rise in concentration during early gestational growth make it an excellent marker for the early detection of pregnancy.

### PRINCIPLE

The test is a one-step, qualitative, double antibodies sandwich lateral flow immunochromatographic assay. The test utilizes a combination of antibodies including a monoclonal hCG antibody to selectively detect elevated levels of hCG. Each test device contains mouse monoclonal anti-β-hCG antibody colloidal gold conjugate pre-dried on a pad. Mouse monoclonal anti-α-hCG antibody (on the Test Line) and goat anti-mouse IgG polyclonal antibody (on the Control Line) are coated and immobilized on a nitrocellulose membrane.

When the absorbent end is immersed into a sample, the sample is absorbed into the device by capillary action and mixes with the antibody-dye conjugate, flowing across the pre-coated membrane.

During the test, hCG in the urine specimen reacts with the dye conjugate (mouse anti-β-hCG antibody- colloidal gold conjugate) and forms a complex. The complex migrates along the membrane to the mouse anti-α-hCG antibody test region (T) and remains captured in the T line.

### PRECAUTIONS

- For in vitro diagnostic use (not for internal use).
- Do not use this test past the indicated expiry date.
- The test device must remain in the sealed pouch until use.
- Do not use if the pouch is torn or not well sealed.
- Do not touch the membrane located within the Result Window and the Absorbent Tip.
- The test should not be reused. The used test should be discarded according to local regulations.
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen

obtained.

- The procedures should be followed precisely for accurate results.

- Keep out of reach of children.

### MATERIALS

#### REAGENTS AND MATERIALS SUPPLIED

Each test kit individually sealed in a foil pouch. Each pouch contains:

- One Pregnancy Test Device
- a desiccant

**NOTE:** The desiccant is only for storage and should not be used during the procedure.

#### MATERIALS REQUIRED BUT NOT SUPPLIED

- Timer,
- A clean container for a urine sample collection (optional)

### STORAGE AND STABILITY

- The test kit is to be stored at 4°C - 30°C in the sealed pouch. **DO NOT FREEZE.** The test is stable through the expiration date printed on the pouch label or label of the outer box.
- Preferably open the pouch only shortly before the test.

### SPECIMEN COLLECTION AND PREPARATION

#### When to do this test?

You can use the pregnancy test from the first day that you have not had your usual period, so from the first day of late menstruation. Note that it is better to wait a few days after the first day of absence of the period, this increases the accuracy of the pregnancy test. There are also so-called "early" tests that can be done a few days before the presumed start date of the period.

It should be known that the earlier the test is performed, the lower the reliability. In addition, it is recommended to conduct a repeat test 5-7 days after the first test to confirm the result.

#### At what time of the day is it better to test?

- Specimens collected at any time of day may be used; however, the first morning urine generally contains the highest concentration of hCG and is recommended for early detection of pregnancy.
- If you do not want to test directly with the urine stream, use a clean, dry glass or plastic container to collect urine.
- Fresh urine does not require any special pre-treatment.

#### Specimen Storage

If the specimen is not to be tested immediately, it may be kept at room temperature for 8 hours or at 2°C to 8°C for up to 72 hours. For prolonged storage, specimens may be frozen and stored at -20°C or below. Do not refreeze. Frozen specimens should be fully thawed and mixed well before testing.

### TEST PROCEDURE

If stored refrigerated, ensure that the pouched test device and specimen are brought to room temperature (about 20 minutes) before use.

- Remove the test device from the sealed pouch and use it as soon as possible.
- Hold the stick by the thumb grip with the exposed


Absorbent Tip pointing downward. Urinate on the Absorbent Tip only until it is thoroughly wet (about 2-3 seconds).

**NOTE:** Do not urinate on the Result Window.

If you prefer, you can urinate into a clean and dry container, then dip only the Absorbent Tip of the midstream test into the urine for at least 10 seconds. Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.

- Lay the test device down on a flat surface.
- Start the time and wait for the colored line(s) to appear. Read the result at 5 minutes. Do not read the results after 15 minutes.

**NOTE:** Depending on the concentration of hCG in the test specimen, positive results may be observed as soon as 10-30 seconds. However, to confirm negative results, the complete reaction time of 5 minutes is required.

### INTERPRETATION OF RESULTS

Pay attention to the position of the Control Line(C) and Test Line (T).

- Positive (+) Result:** Two colored lines appear. One line should be in the Control Line region(C) and another line should be in the Test Line region(T).

**NOTE:** This means that your urine contains hCG and you are probably pregnant. The intensity of the color of the bands may vary depending on the different levels of hCG concentration at various stages of pregnancy. You should consult your doctor. He can confirm your pregnancy after analyzing the laboratory results. A test performed under special conditions may also result in a wrong result. The doctor will be able to give you more information and advise you on your pregnancy. See also " LIMITATIONS ".

- Negative (-) Result:** One colored line in the Control Line region(C). No apparent colored line appears in the Test Line region(T).

**NOTE:** It means that no hCG has been detected in your urine. You may not be pregnant, or the level of pregnancy hormone may not be high enough yet to be detected. If you still have reason to believe that you are pregnant (no period, for example), you can repeat the test a few days later. If the test is still negative, consult your doctor. See also "LIMITATIONS ".

- Invalid Result:** At 5 minutes, if no lines appear, or a test line appears without a Control line, the result is invalid, and the test should be repeated using a new device.

**NOTE:** The intensity of the red color in the Test Line region will vary depending on the concentration of hCG present in the specimen. Therefore, any shade of color in the Test

Line region (T) should be considered positive.

### QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the Control Line region (C) is considered an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative procedural control. If a background color appears in the result window and interferes with the ability to read the test result, the result may be invalid.

### LIMITATIONS

- The contents of this kit are for use in the qualitative detection of hCG in urine only, neither the quantitative value nor the rate of increase in hCG can be determined by this test.
- A false negative result (negative when pregnancy exists) may occur if the urine is too dilute or with a very early stage pregnancy. If pregnancy is still suspected, retest using a fresh urine specimen at least 48 hours after initial testing.
- This test may produce false positive results. The presence of hCG in urine should not be used to diagnose pregnancy unless these conditions have been ruled out.

- Fertility drugs containing hCG can give misleading results (these fertility drugs are usually given by injection and testing too soon after administration may give a false 'Pregnant' result).

- hCG may remain detectable for a few days to several weeks after delivery, spontaneous abortion or termination.
- Several medical conditions other pregnancy, including Ectopic pregnancies, ovarian cysts, menopause can cause elevated levels of HCG, which can result in a false or irregular result.

Several conditions other than pregnancy including trophoblastic disease and certain non-trophoblastic neoplasms, including testicular tumours, prostate cancer, breast cancer and lung cancer, cause elevated levels of hCG.

- A conception in the uterus that is lost in the early days of implantation, such as a chemical pregnancy, might initially give a positive test. When the pregnancy becomes nonviable, usually because of the failure of implantation, the low level of HCG might still be detectable on a blood test for a short while. A repeat pregnancy test will be negative a few days later.

- Alcohol may interfere the test result. It is not recommended using the test after drinking.
- Oral contraceptives, painkillers, antibiotics or hormone therapies that do not contain HCG should not affect the test results.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.

### EXPECTED VALUES

- Negative results are expected in healthy non-pregnant women.
- The amount of hCG will vary greatly with gestational age and between individuals.

### PERFORMANCE CHARACTERISTICS

#### 1. Sensitivity

This sensitivity level has been confirmed with hCG standards in urine. The results demonstrated that the analytical sensitivity of the test (the lowest concentration that yields 100% positive results) is 25 IU/L.

#### 2. Specificity

The test was determined from cross reaction studies with known amounts of luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH). The α subunit of TSH, LH, FSH is like that of hCG, which may cause cross reactivity between those hormones. The addition of LH (500 IU/L), FSH (1,000 IU/L), and TSH (1,000 mIU/L) to hCG negative and positive specimens showed no cross-reactivity.

#### 3. Interference study

The following potentially interfering substances were added to hCG negative and positive specimens.

Substance	Concentration	Substance	Concentration
Glucose	2g/dL	Acetaminophen	20 mg/dL
Albumin	2g/dL	Atropine	20 mg/dL
protein	2g/dL	Aspirin	20 mg/dL
Bilirubin	2mg/dL	Ascorbic Acid	20 mg/dL
Hemoglobin	1mg/dL	Ampicillin	20mg/dL
Vitamins C	20 mg/dL	Salicylic Acid	20 mg/dL
Caffeine	20 mg/dL	Phenothiazine	20 mg/dL
Gentisic Acid	20 mg/dL	Thiophene	20 mg/dL
Tetracycline	20mg/dl	Ephedrine	20mg/dL
Acetylsalicylic Acid	20 mg/dL	Ethanol	1%
hCG β-core fragment	1,000,000 pmol/L	Phenylpropanolamine urine specific gravity	20 mg/dL 1.000~1.035
Urine pH value	3~10		

At the levels tested, none of these substances interfered with the expected test results.

#### 4. Accuracy

A multi-center clinical evaluation was conducted comparing the results obtained using the LuxmedIQ One Step Pregnancy Urine Test to another commercially available urine membrane hCG test (predicate device). The study included 600 urine specimens. The results demonstrated a 100% overall agreement (for an accuracy of > 99%) of the LuxmedIQ One step Pregnancy Urine Test when compared to the other urine membrane hCG test.

#### 5. High dose effect

Normal urine that were spiked with hCG concentrations of 62,500, 125,000, 250,000, 500,000, 1,000,000, and 2,000,000 IU/L were used to study the high dose hook effect on Pregnancy Test. It was noticed that both color bands at the test band region and the control region were visible.

#### 6. Standardization

The test has been standardized to the WHO International Standard 5th WHO Chorionic Gonadotrophin (NIBSC code: 07/364) .

### REFERENCES

- Johnson SR, Miro F, Barrett S, Ellis JE. Levels of urinary human chorionic gonadotrophin (hCG) following conception and variability of menstrual cycle length in a cohort of women attempting to conceive. *Curr Med Res Opin.* 2009 Mar; 25(3):741-8
- Pike J, Ellis J, Johnson S, Shaw R, Parkinson P, Perry P. Accuracy of new Clearblue Digital Pregnancy Test with Conception Indicator. *Hum Reprod* 2010; 25(suppl 1):i169-170

- Cole LA, Ladner DG, Byrn FW. The normal variabilities of the menstrual cycle. *Fertil Steril.* 2008 Apr 21.

### INDEX OF SYMBOLS

	Do not re-use	 35 °C	Temperature Limit
	Catalogue Number		Consult instructions
	In Vitro Diagnostics		European Representative
	Use-by-date		Batch/Lot Code
	Keep dry		Quantity
	Keep away from sunlight		CE Symbol
	Manufacturer		

### CONTACT

**EQUINOX BIOTECH CO., LTD.**  
R116, 1/F, Transportation Bureau Building, No. 95, Yingbin Avenue, Huacheng Subdistrict, Huadu District, Guangzhou 510806, Guangdong, China

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Substance	Concentration	Substance	Concentration
Glucose	2g/dL	Acetaminophen	20 mg/dL
Albumin	2g/dL	Atropine	20 mg/dL
protein	2g/dL	Aspirin	20 mg/dL
Bilirubin	2mg/dL	Ascorbic Acid	20 mg/dL
Hemoglobin	1mg/dL	Ampicillin	20mg/dL
Vitamins C	20 mg/dL	Salicylic Acid	20 mg/dL
Caffeine	20 mg/dL	Phenothiazine	20 mg/dL
Gentisate	20 mg/dL	Thiophene	20 mg/dL
Genistic Acid	20 mg/dL	Ephedrine	20mg/dL
Tetracycline	20mg/dl	Ethanol	1%
Acetylsalicylic Acid	20 mg/dL	Phenylpropanolamine	20 mg/dL
hCG β-core fragment	1,000,000 pmol/L	urine specific gravity	1.000~1.035
Urine PH value	3~10		

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### REFERENZEN

- Johnson SR, Miro F, Barrett S, Ellis JE. Levels of urinary human chorionic gonadotrophin (hCG) following conception and variability of menstrual cycle length in a cohort of women attempting to conceive. *Curr Med Res Opin.* 2009 Mar; 25(3):741-8
- Pike J, Ellis J, Johnson S, Shaw R, Parkinson P, Perry P. Accuracy of new Clearblue Digital Pregnancy Test with Conception Indicator. *Hum Reprod* 2010; 25(suppl 1):i169-170
- Cole LA, Ladner DG, Byrn FW. The normal variabilities of the menstrual cycle. *Fertil Steril.* 2008 Apr 21.

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### KONTAKT

	<b>EQUINOX BIOTECH CO., LTD.</b> R116, 1/F, Transportation Bureau Building, No. 95, Yingbin Avenue, Huacheng Subdistrict, Huadu District, Guangzhou 510806, Guangdong, China		<b>MedNet EC-REP GmbH</b> Borkstrasse 10 48163 Münster, Germany
--	--	--	---

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